



## PRODUCT FEATURES

### SPIROMETRY FILTER 59

For use during pulmonary lung function testing with **CareFusion and Jaeger spirometers**

**Description:** Barrier filtration for protection of patient, staff and equipment from bacterial and viral contamination.

**Code items:**

- individually packaged: **30816**

**Materials:**

- housing: Polystyrene (HIPS) recyclable  
 - filtration membrane: Hydrophobic Polypropylene (SMS)

**Kind of filtration:** Electrostatic / Mechanical

**Flow resistance:** (complies with ATS)

- flow rate 12,00 l/s 0,93 cm H<sub>2</sub>O / (l/s)

**Cross Contamination Reduction Efficiency:** > 99,99 %

**Dead space:** 32,00 ml  
 (connector not included)

**Filtration area:** 25,62 cm<sup>2</sup>  
**Exchange filtration area:** > 72,00 cm<sup>2</sup>

**Weight:** 18 g

**Working procedure:** **CE** 0373 Class IIa - Italian Superior Institute of Health  
 Disposable

**Quantity / box:** 100 pcs.

NOTES: single patient use. PVC and Latex free.

Medical Devices Regulations - Directive 93/42/EEC

**Any size of adaptors is available upon request.** Coronati Consulting: BIOCOMPATIBILITY TESTS



Spirometry Filters